



Cochlear™ Nucleus® 8
Sound Processor

Cochlear™ Nucleus® Kanso® 2
Sound Processor


Hearing with a cochlear implant

How hearing works with a cochlear implant

The Nucleus[®] cochlear implant bypasses parts of the ear that no longer work properly by sending signals directly to the hearing nerve.

- 1 Microphones on the sound processor pick up sounds and the processor converts them into digital information.
- 2 This information is transferred through the coil to the implant just under the skin.
- 3 The implant sends electrical signals down the electrode into the cochlea.
- 4 The hearing nerve fibres in the cochlea pick up the signals and send them to the brain, giving the sensation of sound.
- 5 **Optional** In case there is residual hearing left in the implanted ear, options are available for an additional acoustic component that is worn in the ear canal.

For many, use of a hearing aid in one ear and a cochlear implant in the other can provide enhanced hearing performance.¹ This combination is referred to as bimodal hearing.



How natural hearing works

Sound is perceived naturally by way of air and bone conduction.

- A Sound waves travel through the ear canal and strike the eardrum.
- B These sound waves cause the eardrum and the three bones within the middle ear to vibrate.
- C These vibrations are transferred to the fluids in the inner ear – known as the cochlea – and cause the tiny hair cells in the cochlea to move.
- D The movement of the hair cells produces neural impulses which are sent along the hearing nerve to the brain, where they are interpreted as sound.